# nëo™ aEEG monitor

### simplified Cerebral Function Monitoring (CFM)





The first aEEG monitoring solution optimized for every care-taker in the NICU.

## nëo<sup>™</sup> aEEG monitor

## Simplified Cerebral Function Monitoring (CFM)

The nëo<sup>™</sup> monitor is the easiest-to-use neonatal CFM solution available. Designed specifically to account for NICU workflows - from electrode application, to initiating and reviewing up to 8 channels of aEEG/EEG bedside monitoring with automatic annotation of seizure activity. With 24/7 remote view, compatibility with electronic medical records, and reports, the neo CFM provides clinicians important information about neonatal brain status and easily integrates into the busy NICU's neonatal brain monitoring protocols.



#### Features

- Simple setup and intuitive navigation
- Largest versatility in different electrodes, including exclusive access to shielded headcaps
- Research grade amplifier provides high quality data with less artifacts and less misinterpretation
- Clear markers of suspected seizure activity
- Remote view anytime anywhere
- Integration with patient electronic health records and report generation for efficient patient documentation



**Research Grade Amplifier** provides high input impedance for improved signal quality



Access To Single Lead Electrodes allows compatibility to gold cup, hydrogel, and subcutaneous needle electrodes

#### neonatal **wave**guard<sup>™</sup> caps



	Sizes	Head circumference		Head circumference		Fabric color scheme		
		Min. (cm)	Max. (cm)	Min. (in)	Max. (in)			
С	Child	43	47	16.9	18.5	•	•	•
I	Infant	39	43	15.4	16.9	•	•	•
В	Baby	36	39	14.2	15.4	•	•	•
N5	Neonatal	33	36	13	14.2	•	•	•
N4	Neonatal	31	33	12.2	13	•	•	•
N3	Neonatal	29	31	11.4	12.2	•	•	•
N2	Neonatal	27	29	10.6	11.4	•	•	•
N1	Neonatal	25	27	9.8	10.6	•	•	•

Size indication serves as a guideline only



#### **Technical Specifications**

Touch Screen Monitor					
Weight	4.5 kg		Ор		
Screen size	15", 16:9 Ratio		Мо		
Dimensions	385 x 290 x 45 mm				
Resolution	1920 x 1080 Pixel Full-HD				
Mount	Integrated VESA-100 interface				

Data Acquisition

Actively shielded inputs

86 x 100 x 16 mm 150 - 1000 mV

6

Max. 8

512 Hz

24 bit

>1 GΩ

Software					
Operation System	Windows 10/11 (64bit)				
Monitoring features	Real-time EEG aEEG (computed)				
	Continuous Burst-Suppression-Ratio				
	(BSR)				
	Continuous Inter-Burst-Interval (IBI)				
	Impedance measurement				
	Event markers				
	Online review mode				
	Automatic annotation of seizure activity				
Hardware					
Processor	Intel Core™ i5				
Working Memory	8 GB				
Connectivity	1x USB 3.0				
Storage	mSATA 480GB SSD				
Graphics	Intel® HD Graphics GT-Series				

Power supply Unit Integrated power supply (100-240V)

Coffuerra

### Cybersecurity

**Bipolar** channels

Resolution

Shielding

Amplifier size

Referential channels

Max. Sampling Rate

Input impedance

Input signal range

The **n**ëo<sup>™</sup> Monitor incorporates current best practices to address Cybersecurity hardening. The system's security measures are designed to ensure the safety of the product and to safeguard sensitive information such as patient data.

- 1. Encryption of the hard drive
- 2. Encryption of patient data in the database
- 3. Credential management at key access points

 $\mathbf{n}$ eo<sup>TM</sup> monitor is CE marked as a medical device in the EU, according to MDD 93/42/EEC, CE class IIa, has FDA clearance under 510(k) in the USA, Medical Device License (MDL) issued by Health Canada and is registered in the ARTG in compliance with the Australian TG(MD)R.

 $\ensuremath{^*\text{The}}$  seizure module is not available in all markets. Check with your distributor.

**wave**guard<sup>TM</sup>original and accessories are CE marked medical devices, according to MDR (EU) 2017/745, CE class I, are registered in the ARTG in compliance with the Australian TG(MD)R, have FDA clearance under 510(k) in the USA and Medical Device License (MDL) issued by Health Canada.

Manufactured by eemagine Medical Imaging Solutions GmbH, Berlin, Germany, ISO 13485 certified. ANT Neuro and eemagine are part of the neuromotion group.

For more information about  $\mathbf{n}$ eo<sup>TM</sup> monitor and the regulatory status in your country, contact us at <u>sales@ant-neuro.com</u>.

Your local distributor:



**ANT Neuro b.v.**, Hengelo, The Netherlands, tel: +31 (0) 850 498 175, fax: +31 (0) 850 493 919, internet: www.ant-neuro.com, e-mail: info@ant-neuro.com

Information in this document is subject to change. **www.ant-neuro.com**