



Safe and hygienic swab management

Designed to eliminate risk

Eliminate the unnecessary handling of contaminated swabs with this unique and modernised swab management system. Swabsafe™ allows clinicians to visually keep track of all used consumables during operating times, whilst minimising the re-handling of healthcare waste.

Swab counting is made easy with simply divided container sections and a fully fitted Swabsafe™ lid - which can be secured post-surgery.



**CENTRAL
MEDICAL
SUPPLIES**

Swabsafe™ containers are
produced entirely from recycled
post-consumer plastic (rPET)



Recommended Practice for Management of Surgical Swabs (AfPP 2022)

"As part of an effective risk management strategy, it is recommended that a suitable contained disposal system is in place. If appropriate, the system used should be sealed and the number of contents recorded on the outside of the disposal system."

"At all times during a surgical procedure the perioperative practitioner must be aware of the location of all swabs."

"Used swabs should be accessible to view by the operating surgeon, anaesthetist, and perioperative practitioner where applicable throughout a clinically invasive procedure."

Accountable Items, Swab, Instrument and Needle Count, Standards & Recommendations for Safe Perioperative Practice – Fifth Edition (AfPP 2022)

Ordering Information

Code	Description	QTY
SWAB-910	Small Swabsafe Container	200
SWAB-920	Medium Swabsafe Container	150
SWAB-930	Large Swabsafe Container	75
SWAB-854	Stand Base Only (add required top)	1
SWAB-830	Stand Top - A 2 x Large and 2 x Small/Medium	1
SWAB-831	Stand Top - B 1 x Large and 1 x Small/Medium	1
SWAB-832	Stand Top - C 4 x Small/Medium	1
SWAB-833	Stand Top - D 2 x Large	1
SWAB-834	Stand Top - E 2 x Small/Medium	1
SWAB-820	Swabsafe (SECA) Organ Scales	1
SWAB-810	Swabsafe (Adam) Counting Scales	1

Containers also available via NHS SC:

- FSW411 – Small Swabsafe Container
- FSW412 – Medium Swabsafe Container
- FSW413 – Large Swabsafe Container

Discover more online



Benefits

- Enhances accurate swab counting
- Eliminates the handling of dirty swabs
- Assists swab weighing procedure
- Includes swab-tie receptacle
- Improves patient safety
- Simple and easy to use
- Safe and hygienic



01538 399 541
sales@centralmedical.co.uk
centralmedical.co.uk

